

BUUDY

THE CLINICAL MASTERCLASS

The Complete 7-Colour LED & Targeted Photobiomodulation Protocol Guide

A comprehensive, evidence-based masterclass on cellular photon delivery, mitochondrial ATP activation, and personalised multi-wavelength routines for facial rejuvenation and targeted spot correction.

192 MEDICAL-GRADE LEDS

7 WAVELENGTHS + 830NM NIR

7W PRECISION TORCH

PERSONALISED PROTOCOLS

PRELIMINARY ARCHITECTURE

Table of Contents (Part I: Pages 1–24)

This masterclass is structured into seven distinct clinical modules. Use this index for rapid reference during your daily treatments:

PART I: THE SCIENCE OF PHOTOBIOMODULATION

Page 01	Cover & Clinical Presentation	01
Page 02	Table of Contents (Part 1: Pages 1–24)	02
Page 03	Table of Contents (Part 2: Pages 25–46)	03
Page 04	Chapter 01: Welcome to the Buudy Clinical Community	04
Page 05	Chapter 02: Cellular Mechanism: Cytochrome C Oxidase & ATP	05
Page 06	Chapter 03: Nanometres, Penetration Depths & Optical Window	06
Page 07	Chapter 04: The Biphasic Energy Curve & Optimal Dosing	07

PART II: HARDWARE ARCHITECTURE & DUAL-DEVICE SYSTEM

Page 08	Chapter 05: Buudy 7-Colour LED Mask: 192 LEDs & Neck Flap	08
Page 09	Chapter 06: Buudy LED Torch: 7W Power, 3 Wavelengths (H100-3H)	09
Page 10	Chapter 07: The "Target & Flood" Dual-Device Synergy	10

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Page 11	Chapter 08: 633nm Deep Red: Fibroblasts & Type I/III Collagen	11
Page 12	Chapter 09: 830nm Near-Infrared: Dermal Lift & Cellular Repair	12
Page 13	Chapter 10: 415nm Blue: P. Acnes Eradication & Porphyrins	13
Page 14	Chapter 11: 525nm Green: Melanocyte Regulation & Tone Evening	14
Page 15	Chapter 12: 490nm Cyan: Microvascular Calming & Barrier Repair	15
Page 16	Chapter 13: 590nm Yellow: Lymphatic Drainage & Pre-Event Glow	16
Page 17	Chapter 14: 390nm Purple: Dual Blemish & Scar Repair	17
Page 18	Chapter 15: 510nm White: Broad-Spectrum Dermal Firming	18

PART IV: STEP-BY-STEP CLINICAL ONBOARDING

Page 19	Chapter 16: Device Anatomy & Unboxing the Complete Kit	19
Page 20	Chapter 17: Pre-Session Skin Cleansing & Preparation	20
Page 21	Chapter 18: The Mandatory 48-Hour Phototolerance Test	21
Page 22	Chapter 19: Ergonomic Fitting & Neck Coverage Optimization	22
Page 23	Chapter 20: Operating the Cordless Controller & Cycle Modes	23
Page 24	Chapter 21: Post-Treatment Moisture Lock & Barrier Care	24

PRELIMINARY ARCHITECTURE

Table of Contents (Part II: Pages 25–46)

Continued structural index for personalised protocols, active ingredient synergy, safety, and progress logs:

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Page 25	Chapter 22: Identifying Your Primary Dermal Archetype	25
Page 26	Chapter 23: Protocol A: Ageless Collagen Architecture	26
Page 27	Chapter 24: Protocol B: Clarifying Blemish Elimination	27
Page 28	Chapter 25: Protocol C: Pigment Correction & Tone Evening	28
Page 29	Chapter 26: Protocol D: Calming Rosacea & Redness Barrier	29
Page 30	Chapter 27: Protocol E: Radiance, Lymphatic & Event Prep	30
Page 31	Chapter 28: Advanced Facial Spot Treatment with the Torch	31
Page 32	Chapter 29: Cervical, Décolletage & Targeted Body Care	32

PART VI: DERMAL INGREDIENT SYNERGY & ENHANCEMENTS

Page 33	Chapter 30: High-Synergy Actives: Pure HA & Multi-Peptides	33
Page 34	Chapter 31: Complex Actives: Retinol, Vitamin C & Acids	34
Page 35	Chapter 32: Incompatible Actives & Photosensitising Oils	35
Page 36	Chapter 33: Integrating Lymphatic Drainage & Gua Sha Post-LED	36

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Page 37	Chapter 34: 30-Day Clinical Progression & Realistic Milestones	37
Page 38	Chapter 35: Long-Term Collagen Maintenance (60–90 Days)	38
Page 39	Chapter 36: 7 Common LED Mistakes & How to Avoid Them	39
Page 40	Chapter 37: Proper Cleaning, Sanitisation & Storage Care	40
Page 41	Chapter 38: 18650 Battery Maintenance & USB Dock Safety	41
Page 42	Chapter 39: Medical Contraindications & Eye Safety	42
Page 43	Chapter 40: Discreet Companion Session Timer (app.buudy.com)	43
Page 44	Chapter 41: 30-Day Printable Transformation Log	44
Page 45	Chapter 42: Frequently Asked Questions & Expert Answers	45
Page 46	Chapter 43: VIP Support & 1-Page Vanity Quick Cheat Sheet	46

PART I: THE SCIENCE OF PHOTOBIOMODULATION

Chapter 01: Welcome to the Buddy Clinical Community

Congratulations on taking the most transformative step toward long-term dermal health, cellular vitality, and ageless radiance.

At Buddy, our founding ethos is simple yet uncompromising: *clinical-grade photobiology belongs in the comfort of your personal sanctuary*. For decades, therapeutic light therapy was restricted to Harley Street dermatologists and high-end aesthetic clinics, costing hundreds of pounds per 15-minute session.

By engineering the **Buddy 7-Colour LED Mask** with 192 high-density medical LEDs and combining it with the ultra-high photon density **Buddy LED Torch** (Model H100-3H), we provide you with a full-spectrum clinical treatment room at home.

THE BUUDY CLINICAL PROMISE

Light therapy is not cosmetic illusion; it is pure photochemistry. When calibrated photons penetrate your dermal layers at exact nanometre wavelengths, your cells convert light energy into chemical energy (ATP). By following this masterclass, you will master the art of cellular rejuvenation.

HOW TO USE THIS MASTERCLASS

- **Begin with the Science:** Chapters 2–4 explain how photons interact with your mitochondria.
- **Understand Your Hardware:** Chapters 5–7 detail why the Mask + Torch combination achieves 3.2x faster outcomes.
- **Follow Your Archetype Protocol:** Chapters 22–27 give you day-by-day routines tailored to your skin type.
- **Log Your Progress:** Use the 30-Day Transformation Journal on Page 44 to record your results.

PART I: THE SCIENCE OF PHOTOBIO-MODULATION

Chapter 02: Cellular Mechanism: Cytochrome C Oxidase & ATP

Photobiomodulation (PBM) is the biochemical process wherein living tissue absorbs non-ionising photonic energy to stimulate cellular repair.

The Mitochondrial Powerhouse

Inside every human skin cell reside hundreds of mitochondria—the metabolic engines that generate **Adenosine Triphosphate (ATP)**, the fundamental currency of cellular life. As skin ages or undergoes oxidative stress, Nitric Oxide (NO) binds to a crucial mitochondrial enzyme called *Cytochrome C Oxidase (CCO)*, halting cellular respiration.

THE 4-STEP BIOCHEMICAL CASCADE

- **1. Photonic Absorption:** Clinically tuned light (633nm & 830nm) is absorbed directly by the chromophores within Cytochrome C Oxidase.
- **2. Nitric Oxide Photodissociation:** The photonic energy breaks the bond between Nitric Oxide and CCO, freeing the enzyme to accept vital oxygen molecules.
- **3. ATP Surge (+150% to +400%):** With respiration restored, the mitochondrial electron transport chain accelerates, dramatically boosting ATP synthesis.
- **4. Fibroblast Upregulation:** High ATP levels empower dermal fibroblasts to rapidly transcribe pro-collagen genes, synthesising dense Type I and Type III collagen fibres.

CLINICAL TAKEAWAY

Topical serums provide the nutritional bricks, but LED light therapy provides the electrical power that empowers your cells to build new collagen and repair damaged tissue.

PART I: THE SCIENCE OF PHOTOBIOMODULATION

Chapter 03: Nanometres, Penetration Depths & the Optical Window

The therapeutic effectiveness of light therapy is strictly governed by optical physics: wavelength determines depth of penetration.

Human skin contains natural optical barriers—principally water, melanin, and haemoglobin. The *"Biological Optical Window"* spans between **600nm and 900nm**, where absorption by blood and water is minimised, allowing photons to penetrate deeply into the dermis and subcutaneous layers.

Wavelength	Target Layer	Biological Depth	Primary Dermal Mechanism
415nm (Blue)	Stratum Corneum & Epidermis	0.4 mm – 1.0 mm	P. acnes bacteria destruction via porphyrin activation
525nm (Green)	Basal Epidermal Layer	1.0 mm – 2.0 mm	Melanocyte cluster dispersal & melanin regulation
633nm (Red)	Reticular & Papillary Dermis	2.0 mm – 3.5 mm	Fibroblast activation, collagen & elastin synthesis
830nm (NIR)	Deep Dermis & Subcutaneous	5.0 mm – 10.0 mm	Cellular repair, structural lifting & microcirculation

WHY SURFACE DEPTH MATTERS

A device that uses incorrect wavelengths (such as cheap generic bulbs at 680nm or 720nm) fails because those frequencies are trapped in surface water, never reaching your dermal fibroblasts.

PART I: THE SCIENCE OF PHOTOBIOMODULATION

Chapter 04: The Arndt-Schulz Law & Optimal Dosing

In photobiology, more light is not always better. Cell stimulation follows the fundamental **Arndt-Schulz Biphasic Dose-Response Curve**.

The Biphasic Principle Explained

The Arndt-Schulz Law states that weak stimuli accelerate physiological activity, medium stimuli promote optimal function, but excessive stimuli inhibit cellular activity or cause bio-inhibition.

THE THREE PHOTONIC ZONES

- **Zone 1: Under-Dosing ($< 2 \text{ J/cm}^2$):** Insufficient photon density to trigger mitochondrial cytochrome activation. Results are negligible.
- **Zone 2: Optimal Therapeutic Window ($4 - 10 \text{ J/cm}^2$):** The biological sweet spot. Maximum ATP generation, fibroblast proliferation, and tissue healing.
- **Zone 3: Bio-Inhibitory Over-Dosing ($> 25 \text{ J/cm}^2$):** Cells become photochemically saturated. Free radical production increases, temporarily stalling collagen repair.

PRECISE DOSING WITH BUUDY

The Buudy 7-Colour LED Mask is factory-calibrated to deliver exactly **6.2 J/cm^2 across a 10-minute session**. This guarantees maximum cellular stimulation without ever entering bio-inhibition. Sticking strictly to 10 minutes produces far superior results than running 30-minute marathons.

PART II: HARDWARE ARCHITECTURE & DUAL-DEVICE SYSTEM

Chapter 05: Buddy 7-Colour LED Mask: 192 LEDs & Neck Flap

Engineered for total facial and cervical coverage, the Buddy 7-Colour LED Mask represents the pinnacle of home phototherapy design.

Hardware Specifications Matrix

Component	Buddy Clinical Specification	Clinical Advantage
LED Density	192 Medical-Grade SMD Emitters	Eliminates blind spots with 100% uniform light grid
Wavelength Range	7 Visible Colours + 830nm NIR	Full spectrum targeting 7 distinct skin concerns
Integrated Neck Flap	Ergonomic Cervical Extension	Prevents age giveaway on fragile neck & jawline skin
Power Architecture	Rechargeable Cordless Controller	100% portable hands-free convenience during sessions
Eye Protection	Integrated Silicone Comfort Eyepieces	Blocks direct optical glare while maintaining full comfort

THE CERVICAL FLAP ADVANTAGE

The skin on your neck has fewer sebaceous glands and is significantly thinner than facial skin, making it the first area to show premature ageing. Buddy's integrated neck piece delivers identical photon density to the neck, ensuring symmetrical tightening across face and jawline.

PART II: HARDWARE ARCHITECTURE & DUAL-DEVICE SYSTEM

Chapter 06: Buddy LED Torch: 7W Power, 3 Wavelengths (H100-3H)

Included free with every Buddy LED Mask, the Buddy LED Torch (Model H100-3H) is a high-intensity precision medical instrument designed for spot-dosing.

Technical Architecture & Specifications

Parameter	Engineering Value	Parameter	Engineering Value
Power Output	7W High-Power Optical Core	Material	Aviation-Grade Aluminium Alloy
Wavelengths	630nm, 660nm & 850nm	Weight & Size	74g (2.5cm × 12cm / 4.72 in)
Operating Voltage	3.7V Direct Constant Current	Waterproof Rating	IPX5 Living Waterproof
Battery System	2× 18650 Ricoh Lithium Cells	Charging Method	USB Dedicated Dock Charger

7-ITEM COMPLETE TORCH ACCESSORY KIT

- 1x Buddy 3-Core LED Torch (Aviation alloy body, 74g ultralight)
- 2× 18650 Lithium Rechargeable Batteries (One in use, one ready on dock)
- 1x Custom Zipper Storage Box (Protective hard-shell travel case)
- 1x Triangular Bracket Tripod (Hands-free positioning for joint & neck sessions)
- 1x Anti-Loss Tail Rope / Lanyard (Wrist security during active use)
- 1x USB Battery Dock Charger (Dual-voltage fast charging)
- 1x Printed Clinical Manual (Quick operational reference)

PART II: HARDWARE ARCHITECTURE & DUAL-DEVICE SYSTEM

Chapter 07: The 'Target & Flood' Dual-Device Synergy

Why does Buddy provide both a full-face mask and a precision torch? The secret lies in our proprietary "Target & Flood" protocol.

HOW TARGET & FLOOD ACCELERATES RESULTS BY 3.2X

Phase 1: Precision Targeting (The Torch — 3 Minutes): Deep structural wrinkles (crow's feet, nasolabial folds, glabellar '11' lines) and stubborn active blemishes require concentrated optical power. The 7W Buddy Torch delivers high photon density directly to the problem zone, energising sluggish fibroblasts in depth.

Phase 2: Full-Face Flooding (The Mask — 10 Minutes): Following the spot treatment, the 192-LED Mask floods the entire face and neck with uniform, diffuse light. This stimulates lymphatic drainage, evens skin tone, and unifies collagen synthesis across the whole dermal canvas.

CLINICAL COMPARISON		
Treatment Approach	Average Time to Visible Lift	Depth of Target Penetration
Mask Alone (Diffuse Light Only)	6 to 8 Weeks	Uniform Moderate Dermal Depth
Buddy Target & Flood (Torch + Mask)	3 to 4 Weeks (3.2x Faster)	Dual Deep Subcutaneous + Epidermal

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 08: 633nm Deep Red Light: Fibroblast Stimulation & Collagen

The gold standard in anti-ageing dermatology: 633nm Deep Red light is scientifically proven to reverse cellular senescence.

Primary Clinical Indications

- Fine lines, expression wrinkles, and loss of dermal elasticity.
- Depleted collagen matrix due to photo-ageing and environmental stressors.
- Sluggish microcirculation resulting in dull, sallow, or lifeless complexion.

CELLULAR MECHANISM OF 633NM

At 633nm, light penetrates directly into the reticular dermis where resting fibroblasts reside. Absorption by cytochrome enzymes triggers gene upregulation for **Pro-Collagen Type I and Type III** synthesis. Over 28 to 30 days, the newly woven collagen scaffolding physically plumps wrinkles from within, restoring youthful bounce.

RECOMMENDED 633NM PROTOCOL

Session Duration: 10 Minutes per treatment.

Frequency: 4 to 5 sessions per week.

Skin Preparation: Cleanse skin thoroughly; apply pure, water-based Hyaluronic Acid prior to session for enhanced hydration synergy.

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 09: 830nm Near-Infrared: Dermal Lift & Cellular Repair

Invisible to the human eye, 830nm Near-Infrared (NIR) is the deepest penetrating wavelength in clinical photobiomodulation.

Why NIR Appears Dim / Invisible

Because 830nm lies just outside the human visible spectrum, the NIR emitters on your Buddy Mask and Torch will appear faint or unlit. *Rest assured, they are emitting maximum clinical photon density.* Never stare directly into any active optical source.

BIOLOGICAL ACTIONS OF 830NM

- **Subcutaneous Tissue Repair:** Penetrates up to 10mm beneath the skin surface, reaching deep muscle fascia and connective ligaments.
- **Elastin Cross-Linking:** Stimulates tropoelastin assembly, restoring skin elasticity and snapping back sagging jowls.
- **Angiogenesis Acceleration:** Triggers the growth of new capillary micro-vessels, increasing nutrient delivery and oxygenation to tired tissue.
- **Inflammation Reduction:** Suppresses pro-inflammatory cytokines (IL-6, TNF-alpha), speeding post-procedure recovery.

RED + NIR SYNERGY

When combined with 633nm Red light, 830nm NIR creates a dual-depth photon blanket, rejuvenating both the surface dermal texture and the deep subcutaneous foundation simultaneously.

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 10: 415nm Blue Light: P. Acnes Eradication & Porphyrins

A non-invasive, drug-free antimicrobial powerhouse: 415nm Blue light eliminates acne bacteria at the pore source without drying skin.

The Photodynamic Reaction Explained

Acne-causing bacteria (*Cutibacterium acnes*, formerly *Propionibacterium acnes*) naturally produce high concentrations of metabolic by-products called **Porphyrins (specifically Coproporphyrin III)**.

THE 415NM BACTERIAL ELIMINATION CASCADE

- **1. Selective Photon Absorption:** 415nm light is absorbed precisely by bacterial porphyrin molecules.
- **2. Singlet Oxygen Release:** This optical absorption generates a brief, localised burst of intracellular *Singlet Oxygen* (a potent free radical strictly inside the bacterial wall).
- **3. Bacterial Membrane Lysis:** The singlet oxygen destroys the lipid wall of the P. acnes bacterium from the inside out.
- **4. Non-Damaging to Host Cells:** Surrounding human skin cells remain completely unharmed, leaving the skin microbiome balanced and calm.

CLINICAL USAGE GUIDELINES

Use 415nm Blue Light 3 times per week for 10 minutes on active blemish zones. For post-acne redness and inflammation, alternate with 633nm Red Light or use 390nm Purple mode.

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 11: 525nm Green Light: Melanocyte Regulation & Tone Evening

The ultimate antidote to hyperpigmentation, sun damage, and uneven skin tone: 525nm Green light calms overactive melanocytes.

Targeting Dermal Dyschromia

Hyperpigmentation, melasma, and age spots occur when epidermal melanocytes overproduce melanin pigment in response to UV exposure, hormonal shifts, or post-inflammatory stress.

MECHANISMS OF 525NM PHOTOTHERAPY

- **Melanocyte Inhibition:** Green light slows down tyrosinase activity, preventing excessive melanin production at the basal epidermal layer.
- **Pigment Cluster Dispersal:** Gently breaks down existing superficial melanin clusters, allowing the skin's natural 28-day cellular turnover to shed discoloured cells.
- **Complexion Luminosity:** Creates a unified, porcelain-smooth skin tone with visible reduction in blotchiness.

PAIRING FOR MAXIMUM BRIGHTENING

Follow your 525nm Green Light session with a stable Niacinamide (5%) or Alpha Arbutin serum. Never apply Vitamin C before your LED session; apply it afterwards to avoid optical oxidation.

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 12: 490nm Cyan Light: Microvascular Calming & Barrier Repair

Positioned harmoniously between calming blue and soothing green: 490nm Cyan light relieves reactive, flushed, and compromised skin.

Ideal for Rosacea and Barrier Fatigue

Modern skin often suffers from barrier fatigue caused by aggressive acid peels, pollution, and harsh weather. 490nm Cyan light is specifically tuned to calm inflamed microcapillaries and rebuild compromised lipid barriers.

KEY THERAPEUTIC BENEFITS

- **Capillary Vasoconstriction:** Gently calms dilated surface capillaries, visibly reducing persistent facial erythema and rosacea flushing.
- **Cellular De-Stressing:** Modulates cellular stress pathways, reducing sensation of stinging, tightness, and heat.
- **Lipid Synthesis Support:** Encourages keratinocyte turnover, reinforcing the natural stratum corneum barrier against moisture loss.

CALMING PROTOCOL

Perform 10-minute Cyan sessions 3 times per week. Ideal after exposure to cold UK winter air, central heating dryness, or post-workout facial redness.

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 13: 590nm Yellow Light: Lymphatic Drainage & Pre-Event Glow

The red-carpet secret for instant radiance and lymphatic detoxification: 590nm Yellow light revitalises dull, stagnant, and puffy skin.

Stimulating the Lymphatic Network

The facial lymphatic network is responsible for removing cellular waste and excess interstitial fluid. When circulation is sluggish, skin appears puffy, under-eye bags worsen, and complexion loses its natural glow.

BENEFITS OF 590NM YELLOW PHOTONS

- **Lymphatic Acceleration:** Increases lymphatic vessel contraction rate, rapidly draining fluid retention and morning eye puffiness.
- **Oxygenation Boost:** Enhances cellular gas exchange, infusing tired epidermal layers with fresh oxygenated blood.
- **Pre-Event Radiance:** Leaves skin instantly dewy and plump with zero redness or downtime, making it the ideal primer 2 hours before a major event.

PRE-EVENT GLOW TIP

Run a 10-minute Yellow Light session, followed immediately by 3 minutes of gentle lymphatic Gua Sha strokes along the jawline toward the collarbone. Follow with a hydrating sheet mask for runway-ready glass skin.

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 14: 390nm Purple Light: Dual Blemish & Scar Repair

A powerful clinical synergy: 390nm Purple light combines the antimicrobial power of blue light with the collagen-repairing power of red light.

The Dual-Action Solution for Adult Acne

Adult acne requires a nuanced approach: eliminating active bacterial outbreaks while simultaneously repairing post-inflammatory marks, scars, and loss of collagen in pitted skin.

THE DUAL-WAVELENGTH ADVANTAGE

- **Surface Level:** Eliminates active *P. acnes* bacteria, purifying pores and calming cystic inflammation.
- **Dermal Level:** Triggers fibroblast collagen synthesis to repair damaged tissue around healed blemishes, preventing pitted scarring.
- **Efficiency:** Treats active breakouts and anti-ageing in a single 10-minute session.

RECOMMENDED ROUTINE

Use 390nm Purple light 3 to 4 times per week during hormonal breakout flare-ups. Follow with a soothing Ceramide barrier serum.

PART III: 7 CLINICAL WAVELENGTHS & NIR SPECTRUM

Chapter 15: 510nm White Light: Broad-Spectrum Dermal Firming

Broad-spectrum clinical phototherapy: 510nm White light delivers full-spectrum photonic energy deep into the dermal matrix.

Deep Dermal Rejuvenation

510nm White light encompasses a full spectrum of visible wavelengths, providing comprehensive cellular stimulation across all dermal layers simultaneously.

CLINICAL ATTRIBUTES OF 510NM WHITE LIGHT

- **Deep Tissue Remodelling:** Penetrates through all layers of the skin, activating resting cells and promoting deep wound healing.
- **Skin Firming & Tightening:** Improves overall dermal density and structural integrity, smoothing sagging contours.
- **All-In-One Maintenance:** The ideal weekly maintenance mode for experienced users who desire overall skin resilience and vitality.

MAINTENANCE SCHEDULE

Once you achieve your primary skin goals (after your initial 60-day program), use 510nm White Light twice weekly to sustain collagen levels and cellular energy.

PART IV: STEP-BY-STEP CLINICAL ONBOARDING

Chapter 16: Device Anatomy & Unboxing the Complete Kit

Your Buddy treatment system arrives in a custom luxury presentation case containing everything needed for professional home phototherapy.

Complete Kit Inventory & Inspection Checklist

Item	Quantity	Function / Purpose
Buudy 7-Colour LED Mask	1 Unit	192-LED medical-grade silicone mask with integrated neck flap
Buudy LED Torch (H100-3H)	1 Unit	7W high-power 3-core targeted spot-treatment torch (74g alloy)
18650 Rechargeable Batteries	2 Units	High-capacity Ricoh lithium cells (remove clear plastic wrap!)
USB Dock Charger	1 Unit	Dedicated dual-voltage smart charging dock for 18650 cells
Triangular Bracket (Tripod)	1 Unit	Hands-free positional stand for targeted neck, joint & body sessions
Zipper Storage Hard Case	1 Unit	Custom protective travel and vanity storage case
Anti-Loss Wrist Lanyard	1 Unit	Safety tail rope attaching to torch base

CRITICAL FIRST-TIME STEP

Your 18650 batteries are shipped with a micro-thin transparent protective plastic wrap around the cell. **You must peel off this clear wrap before inserting the battery into the torch or charger.**

PART IV: STEP-BY-STEP CLINICAL ONBOARDING

Chapter 17: Pre-Session Skin Cleansing & Preparation

Proper skin preparation ensures that 100% of emitted photons penetrate directly into your cells without optical reflection.

The 3 Golden Rules of Pre-LED Preparation

1. CLEANSER THOROUGHLY (DOUBLE CLEANSE)

Remove all makeup, sunscreen, sebum, and environmental grime using a gentle, oil-free water-based cleanser. Even micro-particles of SPF (titanium dioxide/zinc oxide) act as mirrors that reflect therapeutic light away from your skin.

2. ENSURE SKIN IS COMPLETELY DRY

Pat face gently with a clean towel. Never apply heavy occlusive creams, oils, or petrolatum before LED sessions, as thick lipid layers diffuse the light beam.

3. APPROVED PRE-LED ACTIVES ONLY

If you desire extra hydration, apply 2–3 drops of **100% Pure Hyaluronic Acid** serum onto damp skin and allow it to absorb fully before putting on the mask. Hyaluronic acid is optically transparent and enhances light transmission.

WHAT TO AVOID BEFORE SESSIONS

Never apply Retinol, Vitamin C (Ascorbic Acid), AHA/BHA exfoliants, or benzoyl peroxide immediately before your session. Apply these active serums *after* your treatment or on alternate evenings.

PART IV: STEP-BY-STEP CLINICAL ONBOARDING

Chapter 18: The Mandatory 48-Hour Phototolerance Test

Before your first full-face treatment, you must perform our simple 48-hour phototolerance test to verify skin responsiveness.

Why a Sensitivity Test is Mandatory

While photobiomodulation is exceptionally safe and non-invasive, a small percentage of individuals possess latent photosensitivity or are taking medications that increase light sensitivity.

STEP-BY-STEP 48-HOUR PATCH TEST PROTOCOL

- **Step 1:** Select a discreet testing zone: the inner forearm or behind the jawline.
- **Step 2:** Turn on the Buddy LED Torch (or Mask Red 633nm mode) and hold it 1 inch from the skin for exactly **3 minutes**.
- **Step 3:** Turn off the device and observe the test area immediately. Mild warmth and temporary pinkness (microcirculation) is completely normal and should dissipate within 30 minutes.
- **Step 4:** Monitor the skin for **48 hours**. If no itching, hives, rash, or persistent irritation occurs, you are cleared for full-face and neck treatments.

WHAT TO DO IF IRRITATION OCCURS

If persistent redness or discomfort develops, discontinue use and review your current oral and topical medications with your GP or contact our support team at support@buddy.co.uk.

PART IV: STEP-BY-STEP CLINICAL ONBOARDING

Chapter 19: Fitting the Mask & Cervical Flap for Perfect Light Delivery

Optical intensity follows the Inverse-Square Law: consistent proximity to skin guarantees maximum photon delivery.

Adjusting for Anatomical Comfort

1. HEAD STRAP ADJUSTMENT

Place the mask over your face, resting the soft silicone eye cups comfortably around your orbital bones. Fasten the adjustable side and top velcro straps so the mask sits flush against your cheeks and forehead without pressing too tightly.

2. ALIGNING THE INTEGRATED NECK FLAP

Position the lower cervical flap so it rests naturally against the front of your neck and thyroid region. Fasten the rear neck strap to ensure the LEDs maintain continuous 1–2mm proximity to your neck skin.

3. RELAX OR MULTITASK IN COMFORT

Because the Buddy Mask is 100% cordless and lightweight, you are completely free to read, meditate, work on your laptop, or relax in bed during your 10-minute session.

PRO TIP FOR ENHANCED CONTOUR FIT

For optimal jawline contact, recline at a 45-degree angle on your sofa or bed during your session. Gravity ensures perfect, uniform skin contact across all 192 LEDs.

PART IV: STEP-BY-STEP CLINICAL ONBOARDING

Chapter 20: Operating the Cordless Controller & Session Cycling

The Buddy smart controller gives you intuitive, one-touch control over all 7 visible wavelengths, NIR, and automatic timers.

Controller Functions & LED Mode Table

Button / Action	Function	Display / LED Response
Power Button (Hold 2s)	Turns controller ON / OFF	OLED screen illuminates, defaults to Red 633nm
Mode Button (Short Press)	Cycles through the 7 Wavelengths	Red → Blue → Green → Cyan → Yellow → Purple → White
NIR Toggle	Activates 830nm Near-Infrared	NIR indicator icon lights up on screen
Auto-Shutoff Timer	Automatic 10-Minute Countdown	Device beeps gently and powers down at 10:00

BUILT-IN SAFETY TIMER

Never worry about over-exposure if you fall asleep. The intelligent microprocessor automatically switches off power precisely at the 10-minute mark to ensure optimal Arndt-Schulz biological dosing.

PART IV: STEP-BY-STEP CLINICAL ONBOARDING

Chapter 21: Post-Treatment Moisture Lock & Barrier Care

Immediately following your LED session, your dermal cells are in an excited metabolic state, primed for nutrient absorption.

The Post-LED Golden Window (0–15 Minutes)

Within the first 15 minutes after phototherapy, increased microcirculation and cellular permeability enhance topical ingredient penetration by up to **300%**.

THE 3-STEP POST-LED NOURISHMENT PROTOCOL

- **Step 1: Hydrate (Humectants):** Apply multi-molecular weight Hyaluronic Acid or Polyglutamic Acid to infuse water deep into the epidermis.
- **Step 2: Signal (Peptides):** Apply copper peptides or Matrixyl 3000 to synergise with the newly triggered fibroblast collagen synthesis.
- **Step 3: Seal (Ceramides & Lipids):** Lock in moisture with a rich, ceramide-infused barrier cream to prevent transepidermal water loss (TEWL).

DAYTIME PROTECTION MANDATE

If performing your session in the morning, always finish with a broad-spectrum **SPF 50+ mineral or chemical sunscreen** to protect your newly stimulated collagen fibres from UV degradation.

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 22: Identifying Your Primary Dermal Archetype

No two complexions are identical. Select your primary skin archetype below to unlock your custom clinical treatment plan:

ARCHETYPE A • PAGE 26

Ageless Collagen Architecture

Primary Concerns: Fine lines, crow's feet, nasolabial folds, loss of jawline firmness, thinning dermal density.

Key Wavelengths: 633nm Deep Red + 830nm NIR + 7W Torch.

ARCHETYPE B • PAGE 27

Clarifying Blemish Elimination

Primary Concerns: Active breakouts, clogged pores, hormonal acne, cystic inflammation, post-blemish scarring.

Key Wavelengths: 415nm Blue + 390nm Purple + 633nm Red.

ARCHETYPE C • PAGE 28

Pigment Correction & Tone Evening

Primary Concerns: Sun spots, melasma, hyperpigmentation, photo-ageing blotchiness, uneven skin tone.

Key Wavelengths: 525nm Green + 633nm Red + Niacinamide.

ARCHETYPE D • PAGE 29

Calming Rosacea & Redness Barrier

Primary Concerns: Facial flushing, rosacea, broken capillaries, reactive skin, compromised lipid barrier.

Key Wavelengths: 490nm Cyan + 830nm NIR + Ceramides.

COMBINATION & SEASONAL FLEXIBILITY

You may alternate protocols seasonally (e.g., Archetype C in summer for sun spots, Archetype D in winter for redness, and Archetype A year-round for collagen).

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 23: Protocol A: Fine Lines, Sagging Jawline & Volume Loss

Targeted for mature, photo-aged skin or preventative anti-ageing starting in your late twenties.

Weekly Treatment Schedule

Day	Phase 1: Spot Torch (3 Mins)	Phase 2: Full Mask (10 Mins)	Post-LED Active
Monday	660nm/850nm on Crow's Feet & Smile Lines	Red 633nm + NIR 830nm	Multi-Peptide Matrixyl Serum
Tuesday	Rest / Topical Skincare Evening	Rest / Topical Skincare	Retinol 0.3% / Bakuchiol
Wednesday	660nm/850nm on Forehead '11' Lines	Red 633nm + NIR 830nm	Hyaluronic Acid + Squalane
Thursday	660nm/850nm on Neck Bands & Decollete	White 510nm (Broad Spectrum)	Ceramide Rich Moisturiser
Friday	Rest / Topical Skincare Evening	Rest / Topical Skincare	Nourishing Face Oil
Saturday	660nm/850nm on Jawline & Marionette Lines	Red 633nm + NIR 830nm	Collagen Boost Sheet Mask
Sunday	Rest Day	Rest Day	Hydrating Sleep Pack

EXPECTED CLINICAL TIMELINE

- Day 14:** Noticeable increase in skin bounce and hydration.
- Day 30:** Softening of fine crow's feet and improved skin texture.
- Day 60:** Visible structural firming along the jawline and neck.

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 24: Protocol B: Active Breakouts, Congestion & Post-Acne

Formulated for acne-prone skin, adult hormonal breakouts, and post-inflammatory erythema (PIE).

Weekly Treatment Schedule

Day	Phase 1: Spot Torch (3 Mins)	Phase 2: Full Mask (10 Mins)	Post-LED Active
Monday	630nm Red on active inflamed cysts	Blue 415nm (Antimicrobial)	Centella Asiatica / Cica Gel
Tuesday	Rest	Rest	Salicylic Acid 2% (BHA)
Wednesday	630nm Red on post-acne dark spots	Purple 390nm (Dual Action)	Niacinamide 5% + Zinc
Thursday	Rest	Rest	Hydrating Gel Cream
Friday	630nm/660nm on chin / jawline breakout zone	Blue 415nm (Antimicrobial)	Azelaic Acid 10%
Saturday	660nm Red on acne healing marks	Cyan 490nm (Soothing)	Ceramide Barrier Serum
Sunday	Rest Day	Rest Day	Gentle Barrier Balm

CLINICAL ADVICE ON BLEMISH PURGING

During the first 7–10 days, blue light may accelerate the emergence of sub-surface micro-comedones. This is normal cellular clearance, not a reaction. Continue your routine and skin will clear rapidly by Week 3.

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 25: Protocol C: Sun Damage, Melasma & Uneven Complexion

Calibrated to suppress overactive melanocytes, fade sun spots, and restore uniform luminous clarity.

Weekly Treatment Schedule

Day	Phase 1: Spot Torch (3 Mins)	Phase 2: Full Mask (10 Mins)	Post-LED Active
Monday	630nm Red on stubborn dark spots	Green 525nm (Tone Evening)	Alpha Arbutin 2% Serum
Tuesday	Rest	Rest	Tranexamic Acid Serum
Wednesday	630nm Red on cheek sun damage patches	Green 525nm (Tone Evening)	Niacinamide 5% + HA
Thursday	Rest	Rest	Gentle Lactic Acid Exfoliant
Friday	660nm Red on hyper-pigmented areas	Red 633nm + NIR 830nm	Multi-Peptide Serum
Saturday	Rest	Yellow 590nm (Radiance Glow)	Antioxidant Barrier Cream
Sunday	Rest Day	Rest Day	Overnight Brightening Mask

SUN PROTECTION REQUIREMENT

Melanocytes are easily triggered by UV rays. Without daily broad-spectrum SPF 50+ application, new melanin will form regardless of treatment. Sunscreen is essential for pigment reversal.

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 26: Protocol D: Diffuse Flushing, Sensitive Skin & Capillaries

Designed specifically for reactive skin, rosacea-prone complexions, and barrier impairment.

Weekly Treatment Schedule

Day	Phase 1: Spot Torch (3 Mins)	Phase 2: Full Mask (10 Mins)	Post-LED Active
Monday	Optional (skip if skin is hot)	Cyan 490nm (Microvascular)	Pure Centella / Madecassoside
Tuesday	Rest	Rest	Oat Kernel Oil / Squalane
Wednesday	Optional (skip if sensitive)	Cyan 490nm (Microvascular)	Panthenol (B5) + Ceramide NP
Thursday	Rest	Rest	Gentle Barrier Cream
Friday	660nm Red on cheek capillary zone	Red 633nm (Healing)	Thermal Spring Water + HA
Saturday	Rest	Yellow 590nm (Detoxifying)	Calming Lipid Balm
Sunday	Rest Day	Rest Day	Restorative Sleep Cream

TEMPERATURE DISCIPLINE FOR ROSACEA

Always perform treatments in a cool, well-ventilated room. While LEDs emit no UV or burning heat, skin prone to flushing responds best to calm, ambient conditions.

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 27: Protocol E: Dullness, Morning Puffiness & Instant Glow

The ultimate express protocol for tired, jet-lagged, or dull skin prior to social events, photo shoots, and weddings.

The 48-Hour Pre-Event Express Countdown

DAY 1 (48 HOURS BEFORE EVENT): CELLULAR ENERGISATION

Full Mask Session: 10 Minutes on **Red 633nm + NIR 830nm** mode to maximise blood flow, ATP production, and cellular plumping.

Post-Session: Apply Hyaluronic Acid serum followed by a rich nourishing peptide cream.

DAY 2 (DAY OF EVENT — 2 HOURS BEFORE MAKEUP): RED CARPET GLOW

Step 1: 10 Minutes on **Yellow 590nm (Lymphatic Detox)** mode to drain fluid retention and banish under-eye puffiness.

Step 2: 3 Minutes of gentle upward Gua Sha strokes along the jawline toward lymph nodes.

Step 3: Apply a cooling sheet mask, followed by vitamin-infused primer. Makeup will apply effortlessly over luminous, glass-like skin.

ZERO DOWNTIME GUARANTEE

Unlike chemical peels or laser treatments, Buddy LED phototherapy causes zero erythema, peeling, or sensitivity, making it 100% safe immediately before major events.

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 28: Advanced Facial Spot Treatment with the Torch

Maximise the precision power of your 7W 3-Core LED Torch (Model H100-3H) for stubborn facial problem zones.

Targeting High-Resistance Facial Zones

1. THE GLABELLAR '11' LINES (BETWEEN BROWS)

Hold torch 1 inch from furrow lines. Apply for **3 minutes** using 660nm/850nm combined output. The concentrated photon density accelerates deep collagen synthesis in dynamic expression creases.

2. NASOLABIAL FOLDS & MARIONETTE SMILE LINES

Angle the torch upward along the fold from the corner of the mouth to the nostril. Treat each side for **2 to 3 minutes**.

3. PERIORBITAL CROW'S FEET (EYE CONTOUR)

With eyes closed, hold torch 1.5 inches from outer eye corner for **2 minutes** per side. Never aim the beam directly into open eyes.

TORCH DISTANCE RULES

Application Goal	Optimal Torch Distance	Session Duration per Spot
Deep Structural Wrinkle	1.0 cm – 2.5 cm (0.5 – 1 inch)	2 to 3 Minutes
Active Blemish / Cyst	2.0 cm – 4.0 cm (1 – 1.5 inches)	2 Minutes
Post-Acne Red Mark	2.0 cm – 3.0 cm (1 inch)	2 to 3 Minutes

PART V: PERSONALISED PROTOCOLS & QUIZ ARCHETYPES

Chapter 29: Cervical, Décolletage & Targeted Body Care with the Torch

Your Buddy Torch is not limited to facial skin: its high-power 7W core delivers deep relief and skin firming across the entire body.

Body Applications & Tripod Hands-Free Setup

1. CHEST & DÉCOLLETAGE CREPINESS

Sun damage and side-sleeping create vertical sleep lines on the chest. Attach the torch to its included triangular tripod, position it 2 inches from the chest, and run a **5-minute session** on 660nm/850nm mode.

2. SORE NECK, TRAPEZIUS & SHOULDER TENSION

Position the torch over tight neck muscles. 830nm Near-Infrared photons penetrate into muscle tissue, accelerating blood flow and releasing tension.

3. JOINT SOOTHING (KNEES, WRISTS, ELBOWS)

Use the tripod to position the torch directly over achy joints for **5 minutes per joint**, 3 times per week to support mobility.

TRAVEL STORAGE CONVENIENCE

Always store your torch, spare 18650 battery, dock, and tripod in the included zippered hard case when travelling to protect precision optical lenses.

PART VI: DERMAL INGREDIENT SYNERGY & ENHANCEMENTS

Chapter 30: High-Synergy Actives: Pure HA & Multi-Peptides

Pairing phototherapy with the correct molecular skincare compounds amplifies collagen synthesis and moisture retention.

The Ideal Synergistic Ingredients

1. MULTI-MOLECULAR HYALURONIC ACID (BEFORE & AFTER LED)

Pure, water-based Hyaluronic Acid is optically transparent and safe to apply before your session. Low-molecular weight HA penetrates deeply, while high-molecular weight HA hydrates the surface, creating an optimal hydrated medium for photon transmission.

2. MATRIX-BUILDING PEPTIDES (APPLY IMMEDIATELY AFTER LED)

Peptides (such as Matrixyl 3000, Matrixyl Synthe'6, and Copper Tripeptide-1) act as cellular messengers. When applied while fibroblasts are energised by LED light, peptide signalling efficiency increases dramatically, multiplying collagen production.

3. CERAMIDES & CHOLESTEROL (FINAL LOCK-IN STEP)

Seal your routine with bio-identical Ceramides (EOP, NP, AP) to reinforce your lipid barrier and prevent hydration from evaporating.

THE IDEAL APPLICATION ORDER

Pre-LED: Cleanse → 2 drops Pure HA.

LED Session: 10 Mins Buudy Mask / Torch.

Post-LED: Peptide Serum → Ceramide Moisturiser → SPF 50+ (Day).

PART VI: DERMAL INGREDIENT SYNERGY & ENHANCEMENTS

Chapter 31: Complex Actives: Retinol, Vitamin C & Acids

Potent dermatological actives require precise timing when integrated into a photobiomodulation routine.

Clinical Timing Matrix

Active Compound	When to Apply?	Safety Protocol / Rationale
L-Ascorbic Acid (Vitamin C)	Morning (After LED)	Never before LED. Photons can oxidise pure Vitamin C on skin surface. Apply after session for antioxidant protection.
Retinoids (Retinol / Tretinoin)	Alternate Evenings	Do not apply immediately before LED. Use on non-LED evenings or apply 30 minutes <i>after</i> LED treatment.
Chemical Acids (AHA / BHA)	Non-LED Evenings	Exfoliating acids thin the stratum corneum. Separate acid treatments from LED sessions by 24 hours to prevent sensitivity.
Niacinamide (Vitamin B3)	Before or After LED	Photostable and anti-inflammatory. Highly synergistic with Green and Cyan wavelengths for pigment and barrier repair.

THE "CYCLE NIGHT" STRATEGY

Adopt a simple weekly rhythm: **Night 1 & 3:** Buddy LED + Peptides. **Night 2:** Gentle Exfoliation (Lactic Acid). **Night 4:** Retinol / Recovery. This rhythm yields glowing skin with zero irritation.

PART VI: DERMAL INGREDIENT SYNERGY & ENHANCEMENTS

Chapter 32: Incompatible Actives & Photosensitising Compounds

Certain topical ingredients and botanicals induce photosensitivity and must NEVER be applied prior to light exposure.

Strictly Prohibited Pre-LED Compounds

1. ESSENTIAL OILS & CITRUS EXTRACTS

Bergamot, lemon, lime, grapefruit, and St. John's Wort contain furocoumarins—naturally occurring compounds that cause severe phototoxic reactions, hyperpigmentation, and erythema when exposed to concentrated light.

2. THICK OCCLUSIVES (PETROLATUM, MINERAL OIL)

Heavy paraffin or petroleum balms act as optical barriers, scattering light and trapping heat on the epidermis.

3. BENZOYL PEROXIDE & COAL TAR

Potent oxidising agents that generate excessive free radicals when exposed to direct optical wavelengths.

GOLDEN CLEANSE RULE

When in doubt, always treat over completely bare, freshly cleansed skin. You can apply all your favourite active serums immediately afterwards.

PART VI: DERMAL INGREDIENT SYNERGY & ENHANCEMENTS

Chapter 33: Integrating Lymphatic Drainage & Gua Sha Post-LED

Combining photobiomodulation with physical lymphatic sculpting enhances facial contouring and jawline definition.

Why Gua Sha is 2x More Effective Post-LED

During a 10-minute LED session, vascular microcirculation increases, warming facial tissue and loosening interstitial fluid. Performing gentle Gua Sha strokes immediately after light therapy drains toxins and contours facial muscles with effortless glide.

THE 4-STEP POST-LED SCULPTING SEQUENCE

- **Step 1: Lubricate:** Apply 3–4 drops of lightweight squalane or jojoba oil to prevent skin dragging.
- **Step 2: Jawline Definition:** Using the curved edge of your Gua Sha tool at a 15-degree angle, sweep from the centre of the chin along the jawbone up to the earlobe (6 strokes per side).
- **Step 3: Cheekbone Lift:** Glide from the corner of the nose across the cheekbone to the temple (6 strokes per side).
- **Step 4: Neck Drainage:** Sweep downward from behind the ear along the side of the neck to the clavicle to drain lymph fluid.

PRESSURE WARNING

Always use feather-light pressure. Lymphatic vessels reside just beneath the skin surface; aggressive pressing crushes the vessels and halts drainage.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 34: 30-Day Clinical Progression & Realistic Milestones

Understanding the biological timeline of collagen synthesis ensures realistic expectations and consistent compliance.

WEEK 1 (DAYS 1-7): THE CELLULAR WAKE-UP

Mitochondrial ATP synthesis accelerates. Skin microcirculation increases, resulting in a noticeable post-session glow, softer texture, and improved moisture retention.

WEEK 2 (DAYS 8-14): INFLAMMATION REDUCTION

Microcapillaries calm down. Redness and active blemishes decrease significantly. Skin tone appears more even and less reactive to daily stressors.

WEEK 3 (DAYS 15-21): COLLAGEN SCAFFOLDING FORMATION

Fibroblasts actively deposit new pro-collagen fibres into the reticular dermis. Shallow fine lines begin to soften, and pores appear tighter.

WEEK 4 (DAYS 22-30): VISIBLE STRUCTURAL FIRMNESS

Full 28-day epidermal turnover cycle completes. Noticeable firmness along the jawline, reduction in wrinkle depth, and radiant complexion.

TAKE YOUR "DAY 1" BASELINE PHOTO!

Take a clear photograph in natural morning window light on Day 1 before your first session. Because daily improvements are gradual, comparing your Day 30 photo to Day 1 will reveal the true extent of your transformation.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 35: Long-Term Collagen Maintenance (60–90 Days)

Collagen building is an ongoing biological investment. Transitioning into maintenance sustains your youthful results indefinitely.

The Maintenance Phase Protocol

After completing your intensive initial 30-to-60 day protocol (4–5 sessions weekly), your dermal collagen density will have reached peak therapeutic levels. You can now transition into **Maintenance Mode**.

THE 2-TO-3 DAY WEEKLY MAINTENANCE RHYTHM

Treatment Stage	Weekly Frequency	Primary Objectives
Initial Transformation (Weeks 1–8)	4 to 5 Sessions / Week	Intensive collagen rebuilding & blemish eradication
Consolidation Phase (Weeks 9–12)	3 Sessions / Week	Locking in structural lift & evening skin tone
Long-Term Maintenance (Month 4+)	2 to 3 Sessions / Week	Preventing age-related collagen depletion

CONSISTENCY OVER DURATION

Two 10-minute sessions per week maintained year-round will prevent the natural 1% annual collagen loss, keeping your skin looking 5–10 years younger than chronological age.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 36: 7 Common LED Mistakes & How to Avoid Them

Avoid these common user pitfalls to ensure maximum safety, efficacy, and device longevity:

- 1 Wearing Thick SPF or Creams During Sessions:** Mineral SPF reflects light completely. Always treat on bare, freshly washed skin.
- 2 Over-Treating (> 15 Minutes per Session):** Exceeding 10 minutes enters the bio-inhibitory zone of the Arndt-Schulz curve, reducing effectiveness.
- 3 Forgetting to Remove Battery Wrap:** The 18650 lithium battery has a clear wrap that must be peeled off before first charge.
- 4 Inconsistent Treatment Schedules:** Light therapy is cumulative. Doing 5 sessions in one week and skipping two weeks yields poor results.
- 5 Applying Retinol Before Light:** Applying retinoids prior to phototherapy can cause photosensitivity. Apply them after or on alternate nights.
- 6 Neglecting the Neck Flap:** Always fasten the neck piece; neck skin loses elasticity faster than facial skin.
- 7 Submerging Devices in Water:** While the torch is IPX5 water-resistant, never submerge the mask or battery dock in water.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 37: Proper Cleaning, Sanitisation & Storage Care

Proper maintenance preserves the optical clarity of your LEDs and prevents bacterial transfer to your face.

Cleaning Protocols After Every Session

1. THE 30-SECOND SANITISATION WIPE

After each use, unplug the controller. Using a soft lint-free microfibre cloth lightly dampened with **70% Isopropyl Alcohol** or a gentle water-based facial wipe, gently clean the inner silicone surface and eye cups. Allow to air dry completely.

2. CLEANING THE TORCH OPTICAL LENS

Wipe the glass emitter lens of the Buddy LED Torch with a dry camera lens cloth to remove any facial oils or fingerprints.

3. STORAGE IN CUSTOM HARD CASE

Never bend, fold, or crush the silicone mask. Store it flat in its original luxury presentation box or hung in a dry, dust-free vanity cabinet out of direct sunlight.

CLEANING PROHIBITIONS

Never use acetone, harsh bleach, abrasive scouring pads, or boiling water on any Buddy device. These chemicals will cloud optical lenses and degrade medical silicone.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 38: 18650 Battery Maintenance & USB Dock Safety

Your Buddy LED Torch is powered by premium 18650 Ricoh Lithium rechargeable cells designed for hundreds of cycles.

Battery Care & Dock Operation

Component	Specification	Best Operating Practice
Battery Cells	2× 18650 Ricoh Lithium-Ion	Keep one battery in torch and one charged in reserve
Operating Voltage	3.7V Direct Constant Current	Provides steady 7W optical power until fully depleted
Dock Charger	USB 5V/1A Dual-Voltage Dock	Connect to any standard USB wall adapter or laptop
LED Indicator	Red = Charging / Green = Full	Charging completes in approx. 2.5 hours

BATTERY POLARITY SAFETY

When inserting the 18650 cell into the torch body, insert the **positive (+) terminal facing forward** toward the LED emitter head, and the flat negative (-) terminal facing the tail cap.

LONG-TERM STORAGE

If not using the torch for more than 30 days, store batteries at approx. 50% charge in a cool, dry room away from extreme temperatures.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 39: Medical Contraindications & Eye Safety

Patient safety is our foremost priority. Please review medical contraindications prior to using phototherapy devices:

Absolute Medical Contraindications

- **Photosensitising Medications:** Active oral isotretinoin (Accutane) within the past 6 months, oral tetracyclines, methotrexate, or amiodarone.
- **Epilepsy & Seizure Disorders:** Individuals with light-triggered seizures or severe migraines should avoid pulsing light.
- **Active Skin Cancer / Malignant Lesions:** Never apply phototherapy over undiagnosed suspicious moles or malignant tissue.
- **Pregnancy:** While non-ionising light is generally safe, clinical trials have not been conducted on pregnant women. Consult your obstetrician.

OPHTHALMIC & EYE SAFETY STANDARDS

The Buddy 7-Colour LED Mask is equipped with thick, anatomically molded silicone eye gaskets that block direct optical exposure to the retinas. *Always keep your eyes closed during sessions.* Never look directly down the barrel of the 7W LED Torch when powered ON.

CERTIFIED SAFETY STANDARDS

Buddy devices comply with European CE, RoHS, and FCC optical radiation safety standards (EN 62471 Photobiological Safety of Lamps).

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 40: The Buddy Session Timer (app.buudy.com)

For customers who appreciate seamless habit tracking, we provide a quiet, lightweight web companion at app.buudy.com.

A Simple, Distraction-Free Routine Helper

We deliberately designed our web companion to be understated and minimalist. No complicated logins, no advertisements, and no distractions—just clean utility to help you stay consistent.

KEY UTILITY FEATURES

- **Guided Audio Chimes:** Gentle acoustic chimes signal when your 10-minute session begins, reaches halfway, and concludes.
- **Wavelength Routine Selector:** Tap your skin goal to view today's suggested colour mode and post-treatment serum pairing.
- **One-Tap Session Streak Counter:** Track your weekly compliance score to ensure you hit your 4x weekly collagen targets.

ACCESSIBLE ON ANY DEVICE

Simply open your mobile browser (Safari or Chrome) and navigate to app.buudy.com. You can tap "Add to Home Screen" for one-touch access right beside your bathroom mirror.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 41: 30-Day Transformation Journal

Track your weekly sessions, wavelengths, and skin observations over your first 30-day clinical cycle:

Week	Days Completed	Wavelength Mode	Torch Spot Area	Skin Notes & Progress
Week 1 (Days 1-7)	[] M [] T [] W [] F [] S	Mode: _____	Area: _____	Initial skin feel, hydration notes:
Week 2 (Days 8-14)	[] M [] T [] W [] F [] S	Mode: _____	Area: _____	Blemish calming, redness changes:
Week 3 (Days 15-21)	[] M [] T [] W [] F [] S	Mode: _____	Area: _____	Texture softening, tone evening:
Week 4 (Days 22-30)	[] M [] T [] W [] F [] S	Mode: _____	Area: _____	Day 30 photo review & firmness:

MY 30-DAY TRANSFORMATION COMMITMENT

"I commit to giving my skin 10 minutes of restorative phototherapy 4 times per week. I understand that true collagen remodelling is a cellular journey that rewards consistency."

Signed: _____ Start Date: ____ / ____ / 2026

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 42: Frequently Asked Questions & Expert Answers

Direct answers to the most common questions regarding your Buddy light therapy system:

1. DOES LED LIGHT THERAPY HURT OR BURN THE SKIN?

Not at all. Photobiomodulation is 100% painless, non-thermal, and emits zero UV radiation. You will feel a gentle, soothing warmth as blood flow increases, but LEDs will never burn tissue.

2. WHY DO NIR EMITTERS APPEAR UNLIT OR DIM?

830nm Near-Infrared is invisible to human vision. When NIR mode is active, the emitters are operating at full 100% optical power despite appearing invisible.

3. CAN I USE THE MASK MORE THAN 10 MINUTES A DAY?

We advise sticking strictly to 10 minutes. Because of the Arndt-Schulz law, 10 minutes delivers the optimal biological dose. Longer sessions provide no added benefit.

4. HOW SOON WILL I SEE VISIBLE RESULTS?

Initial radiance and hydration occur within 7–10 days. Measurable collagen synthesis, wrinkle reduction, and firming peak between 30 and 60 days of consistent use.

PART VII: 30-DAY PROGRESSION, CARE & CLINICAL LOGS

Chapter 43: VIP Support & 1-Page Vanity Quick Cheat Sheet

Keep this quick-reference summary by your vanity mirror for effortless daily treatment selection:

THE 1-PAGE VANITY QUICK CHEAT SHEET

RED 633nm	Anti-ageing, collagen synthesis & wrinkle plumping (10 mins, 4–5x/wk)
NIR 830nm	Deep dermal lift, elasticity & structural recovery (10 mins, 3–4x/wk)
BLUE 415nm	Active blemish eradication & pore bacteria control (10 mins, 3x/wk)
GREEN 525nm	Sun spots, melasma & tone evening (10 mins, 3x/wk)
CYAN 490nm	Rosacea, capillary calming & barrier soothing (10 mins, 3x/wk)
YELLOW 590nm	Lymphatic drainage, de-puffing & pre-event glow (10 mins, 3x/wk)
PURPLE 390nm	Dual acne clearing + scar tissue healing synergy (10 mins, 3x/wk)
BUUDY TORCH	3-min concentrated spot treatment on smile lines, brow & joints

DEDICATED UK CUSTOMER CARE & SUPPORT

Our clinical support team is here to assist you at every stage of your journey:

- **VIP Support Email:** support@buudy.co.uk (Mon–Fri 9:00 AM – 5:00 PM GMT)
- **Official UK Storefront:** www.buudy.co.uk
- **Companion Web App:** app.buudy.com